

RF Group
Group Procedure RF005
Original Issue Date: 07/30/1998
Revision 00

BROOKHAVEN NATIONAL LABORATORY

Procedure for Repairing _____ RF System

Any time the RF Group is going to perform any type of work repairs, modification to an existing or new RF system that involves water in the AGS or Booster ring the mechanical services group is going to assist the RF group with the water.

Head of Safety: _____

Check Each Step

- | | |
|---------------------------------------|-------|
| 1. Secure Water | _____ |
| 2. Drain Water into proper container | _____ |
| 3. RF repairs | _____ |
| 4. Refill using liquid from container | _____ |
| 5. Check system for leaks | _____ |
| 6. Top off System | _____ |
| 7. Turn System on | _____ |

Mechanical Services Group	_____	_____
	Signature	Date

RF Group	_____	_____
	Signature	Date